

MicroRNA profiling of follicular lymphoma identifies microRNAs related to cell proliferation and tumor response


Weixin Wang,¹ Meghan Corrigan-Cummins,¹ Justin Hudson,² Irina Maric,¹ Olga Simakova,¹ Sattva S. Neelapu,³ Larry W. Kwak,³ John E. Janik,⁴ Barry Gause,⁵ Elaine S. Jaffe,⁶ and Katherine R. Calvo¹

¹Hematology Section, Department of Laboratory Medicine, National Institutes of Health Clinical Center, Bethesda, MD; ²Howard University College of Medicine, Washington, DC; ³Department of Lymphoma and Myeloma, Division of Cancer Medicine, The University of Texas M.D. Anderson Cancer Center, Houston, TX; ⁴Metabolism Branch, National Cancer Institute, Bethesda, MD; ⁵SAIC, National Cancer Institute (NCI), Frederick, MD; and ⁶Laboratory of Pathology, National Cancer Institute, Bethesda, MD, USA

Citation: Wang W, Corrigan-Cummins M, Hudson J, Maric I, Simakova O, Neelapu SS, Kwak LW, Janik JE, Gause B, Jaffe ES, and Calvo KR. MicroRNA profiling of follicular lymphoma identifies microRNAs related to cell proliferation and tumor response. Haematologica 2012;97(4):586-594. doi:10.3324/haematol.2011.048132







Online Supplementary Figure S1. miR-20a, -20b, and -194 expression in OCI-Ly8 cells. The levels of endogenous expression of miR-20a, -20b, and -194 were determined by TaqMan quantitative real-time PCR assay in triplicate experiments. RNA from seven FH B cells was used as the control. Four small nuclear and nucleolar RNA, RNU66, RNU38b, Z30, and RNU6B, were used as the internal control. The Δ Ct value was calculated as Ct value of target miR probe subtracted from the geometrical mean of Ct values of RNU66, RNU38b, Z30, and RNU6B. The relative miR expression level was calculated as $2^{\Delta Cl}$. Mean and standard deviation of triplicate experiments are shown for OCI-Ly8 cells; the line indicates the mean value for seven FH B-cell samples. P=0.006, 2.16-fold for miR-20a; P=0.006, 0.52-fold for miR-20b; and P=0.006, 1.86-fold for miR-194.

| ABCC11 | CXCR5 | IRF4 | PTEN |
|---------|----------|----------|-----------|
| ABCE1 | DGCR8 | ITK | RAD51 |
| ABL1 | DICER1 | JAK1 | RECQL |
| ABL2 | DOCK2 | JAK2 | REL |
| ADAP1 | DROSHA | JUN | REV1 |
| AICDA | E2F1 | KDM6B | ROR1 |
| AKT1 | E2F5 | KIAA1274 | RPS6KB1 |
| API5 | EIF2C1 | KLF12 | RUNX1 |
| ATM | EIF2C2 | KLF4 | S100A9 |
| ATR | EIF2C4 | KLK10 | SERPINA9 |
| B3GAT1 | EIF4EBP1 | LAT | SIRT1 |
| BACH1 | EP300 | LEF1 | SMAD1 |
| BACH2 | EPHA2 | LGALS3 | SMAD4 |
| BAX | ERBB3 | LMO2 | SMAD5 |
| BCL11B | ERCC5 | LTA | SOCS1 |
| BCL2 | EZH2 | MAF | SOCS2 |
| BCL6 | FAS | MALT1 | SOX2 |
| BLM | FASLG | MAPK1 | SPI1 |
| BRCA1 | FCER2 | МАРК3 | SPIB |
| BRCA2 | FOXM1 | MDM2 | STAT1 |
| C4A | FOXP1 | MED1 | STAT2 |
| CASP1 | FOXP3 | MIDN | STAT3 |
| CASP2 | FYN | MKI67 | STAT4 |
| CASP8 | GADD45A | MLL | STAT5A |
| CAV1 | GADD45B | MLL2 | SYK |
| CCL2 | GAPDH | MME | TAB2 |
| CCND1 | GCET2 | MMP9 | TANK |
| CCND2 | GUSB | MTOR | TERT |
| CCNL2 | HDAC1 | MYB | TGFB1 |
| CCR7 | HGF | MYC | TNF |
| CD38 | HOXB2 | NEK2 | TNFRSF13B |
| CD3E | HOXD8 | NEUROG1 | TNFRSF17 |
| CD40 | HPRT1 | NFKB1 | TNFRSF1A |
| CD40LG | IBTK | NFKB2 | TNFRSF1B |
| CD7 | IFNG | NFKBIA | TNFRSF25 |
| CD8B | IKBKB | NOTCH1 | TNFSF10 |
| CDC20 | IKBKG | NUCB1 | TNFSF12 |
| CDC42 | IL10 | PAK4 | TNFSF13 |
| CDK6 | IL10RA | PAX5 | TNFSF13B |
| CDKN1A | IL15 | PDCD1 | TNFSF9 |
| CEBPB | IL1B | PGK1 | TP53 |
| CFLAR | IL21 | PIK3R1 | TRAF1 |
| CHEK1 | IL21R | POU3F2 | TXNIP |
| COL9A1 | IL4 | PPM1E | UBE2I |
| CSNK1G2 | IL4R | PRDM1 | USP14 |
| CTNNB1 | IL6 | PRKCB | VEGFA |
| CUX1 | IL7R | PRKCE | XBP1 |
| CXCL12 | IL8 | PRKCQ | ZEB1 |
| CXCL13 | ING1 | PRKCZ | ZEB2 |
| CXCR4 | INPPL1 | PRKRA | |
| | | | |

| Systematic Name | FL vs FH-regulation | p-value | Fold change |
|-----------------|---------------------|---------|-------------|
| hsa-let-7a | up | 0.02 | 13.55 |
| hsa-let-7c | up | 0.02 | 9.93 |
| hsa-let-7f | up | 0.02 | 18.60 |
| hsa-let-7f-1* | down | 0.02 | -2.00 |
| hsa-miR-106a* | down | 0.02 | -2.87 |
| hsa-miR-1205 | down | 0.02 | -2.28 |
| hsa-miR-1224-3p | down | 0.02 | -2.20 |
| hsa-miR-1243 | down | 0.02 | -2.36 |
| hsa-miR-1251 | down | 0.02 | -2.02 |
| hsa-miR-1252 | down | 0.02 | -2.06 |
| hsa-miR-1260 | up | 0.02 | 7.26 |
| hsa-miR-1274a | up | 0.02 | 11.61 |
| hsa-miR-1274b | up | 0.02 | 9.14 |
| hsa-miR-128 | up | 0.02 | 5.90 |
| hsa-miR-1287 | up | 0.01 | 9.49 |
| hsa-miR-1288 | up | 0.02 | 5.08 |
| hsa-miR-1295 | up | 0.02 | 8.35 |
| hsa-miR-1297 | down | 0.02 | -2.20 |
| hsa-miR-1322 | down | 0.02 | -2.40 |
| hsa-miR-135a | down | 0.02 | -2.14 |
| hsa-miR-138-1* | down | 0.03 | -2.03 |
| hsa-miR-141 | down | 0.01 | -6.26 |
| hsa-miR-141* | down | 0.03 | -2.02 |
| hsa-miR-143* | down | 0.02 | -2.10 |
| hsa-miR-147 | down | 0.02 | -2.23 |
| hsa-miR-1471 | up | 0.01 | 9.00 |
| hsa-miR-151-5p | up | 0.02 | 8.26 |
| hsa-miR-15b | up | 0.02 | 14.30 |
| hsa-miR-16 | up | 0.02 | 7.79 |
| hsa-miR-17* | down | 0.01 | -3.49 |
| hsa-miR-187 | down | 0.02 | -2.11 |
| hsa-miR-188-3p | down | 0.02 | -2.38 |
| hsa-miR-18b* | down | 0.02 | -2.83 |
| hsa-miR-1910 | down | 0.02 | -2.04 |
| hsa-miR-1913 | down | 0.02 | -2.42 |
| hsa-miR-193a-5p | up | 0.01 | 8.79 |
| hsa-miR-193b* | up | 0.01 | 7.66 |
| hsa-miR-194 | up | 0.02 | 5.85 |
| hsa-miR-195 | up | 0.02 | 8.53 |
| hsa-miR-19a* | down | 0.02 | -2.14 |
| hsa-miR-200a | down | 0.03 | -2.87 |
| hsa-miR-200b | down | 0.02 | -5.20 |
| hsa-miR-202 | down | 0.01 | -3.30 |
| hsa-miR-203 | down | 0.02 | -6.80 |
| hsa-miR-205 | down | 0.01 | -31.12 |
| | | | 32.22 |

continued on next page

| hsa-miR-20a | up | 0.02 | 9.09 |
|------------------|------|--------------|-----------|
| hsa-miR-20b | up | 0.02 | 7.27 |
| hsa-miR-218-1* | down | 0.02 | -2.22 |
| hsa-miR-219-1-3p | down | 0.02 | -2.28 |
| hsa-miR-22 | down | 0.02 | -2.06 |
| hsa-miR-221* | up | 0.02 | 6.00 |
| hsa-miR-222 | down | 0.02 | -2.12 |
| hsa-miR-222* | down | 0.02 | -2.39 |
| hsa-miR-26a-2* | down | 0.02 | -2.20 |
| hsa-miR-26b | up | 0.02 | 13.16 |
| hsa-miR-27b* | down | 0.02 | -2.00 |
| hsa-miR-297 | down | 0.02 | -2.20 |
| hsa-miR-299-5p | down | 0.02 | -2.12 |
| hsa-miR-301b | down | 0.02 | -2.58 |
| hsa-miR-302d* | down | 0.02 | -2.07 |
| hsa-miR-30a | down | 0.01 | -3.09 |
| hsa-miR-30b | up | 0.02 | 15.15 |
| hsa-miR-323-5p | down | 0.02 | -2.02 |
| hsa-miR-323b-5p | down | 0.02 | -2.64 |
| hsa-miR-330-5p | down | 0.02 | -2.07 |
| hsa-miR-33a | down | 0.01 | -3.22 |
| hsa-miR-345 | up | 0.01 | 6.78 |
| hsa-miR-34a | up | 0.02 | 17.50 |
| hsa-miR-363* | down | 0.02 | -2.08 |
| hsa-miR-374a | up | 0.02 | 23.72 |
| hsa-miR-374b | up | 0.02 | 15.79 |
| hsa-miR-379* | down | 0.02 | -2.08 |
| hsa-miR-380 | down | 0.02 | -2.08 |
| hsa-miR-410 | down | 0.02 | -2.36 |
| hsa-miR-429 | down | 0.02 | -2.16 |
| hsa-miR-431* | down | 0.02 | -3.27 |
| hsa-miR-449b | down | 0.02 | -2.12 |
| hsa-miR-450b-3p | down | 0.02 | -2.54 |
| hsa-miR-451 | up | 0.02 | 7.52 |
| hsa-miR-454 | up | 0.02 | 9.02 |
| hsa-miR-454* | down | 0.02 | -2.22 |
| hsa-miR-486-3p | up | 0.02 | 5.60 |
| hsa-miR-497* | down | 0.02 | -2.01 |
| hsa-miR-502-3p | up | 0.02 | 6.23 |
| hsa-miR-513a-5p | up | 0.01 | 15.39 |
| hsa-miR-513b | up | 0.01 | 7.86 |
| hsa-miR-520h | down | 0.02 | -2.08 |
| hsa-miR-532-3p | up | 0.02 | 5.94 |
| hsa-miR-544 | down | 0.02 | -2.05 |
| hsa-miR-548b-3p | down | 0.02 | -2.02 |
| hsa-miR-548h | down | 0.02 | -2.36 |
| hsa-miR-548j | down | 0.02 | -2.05 |
| hsa-miR-558 | down | 0.02 | -2.23 |
| | | continued on | nort naga |

continued on next page

continued from previous page

| | ¥ | | |
|---------------------------------|------|---------------------------|-------|
| hsa-miR-562 | down | 0.02 | -2.15 |
| hsa-miR-567 | down | 0.02 | -2.02 |
| hsa-miR-569 | down | 0.02 | -2.15 |
| hsa-miR-570 | down | 0.02 | -3.31 |
| hsa-miR-574-3p | up | 0.01 | 5.96 |
| hsa-miR-582-3p | down | 0.02 | -2.28 |
| hsa-miR-584 | up | 0.01 | 11.69 |
| hsa-miR-591 | down | 0.02 | -2.16 |
| hsa-miR-599 | down | 0.02 | -2.04 |
| hsa-miR-600 | down | 0.02 | -2.57 |
| hsa-miR-603 | down | 0.02 | -2.69 |
| hsa-miR-604 | down | 0.02 | -2.42 |
| hsa-miR-606 | down | 0.02 | -2.44 |
| hsa-miR-614 | down | 0.02 | -2.38 |
| hsa-miR-615-3p | down | 0.02 | -2.27 |
| hsa-miR-626 | down | 0.02 | -2.24 |
| hsa-miR-628-5p | down | 0.02 | -2.01 |
| hsa-miR-641 | down | 0.02 | -2.34 |
| hsa-miR-646 | down | 0.02 | -2.12 |
| hsa-miR-647 | down | 0.02 | -2.25 |
| hsa-miR-653 | down | 0.02 | -2.12 |
| hsa-miR-654-3p | down | 0.02 | -2.38 |
| hsa-miR-655 | down | 0.02 | -2.24 |
| hsa-miR-660 | up | 0.02 | 9.84 |
| hsa-miR-661 | down | 0.02 | -2.03 |
| hsa-miR-663 | up | 0.01 | 9.88 |
| hsa-miR-664 | up | 0.01 | 7.69 |
| hsa-miR-664* | up | 0.01 | 9.65 |
| hsa-miR-708* | down | 0.02 | -2.47 |
| hsa-miR-7-1* | up | 0.02 | 7.09 |
| hsa-miR-767-5p | down | 0.02 | -2.27 |
| hsa-miR-873 | down | 0.02 | -2.17 |
| hsa-miR-888 | down | 0.02 | -2.40 |
| hsa-miR-892a | down | 0.02 | -2.05 |
| hsa-miR-9 | up | 0.02 | 8.71 |
| hsa-miR-9* | up | 0.02 | 10.46 |
| hsa-miR-922 | down | 0.02 | -2.32 |
| hsa-miR-935 | down | 0.02 | -2.22 |
| hsa-miR-96 | up | 0.02 | 8.88 |
| hsa-miR-96* | down | 0.02 | -2.25 |
| hsa-miR-98 | up | 0.02 | 7.64 |
| eresioned, it initial terms (T) | 02°E | , -, -, -, - ; | **** |

| GENE ID | FL vs FH-regulation | p-value | Fold change |
|-------------|---------------------|---------|-------------|
| ADAP1 | down | 0.05 | -4.79 |
| AKT1 | up | 0.02 | 1.72 |
| BACH1 | down | 0.03 | -2.03 |
| BAX | up | 0.03 | 1.82 |
| BCL2 | up | 0.02 | 16.29 |
| BCL6 | up | 0.03 | 1.91 |
| BRCA1 | down | 0.04 | -2.14 |
| CDKN1A | down | 0.03 | -11.80 |
| CHEK1 | down | 0.02 | -4.15 |
| CUX1 | up | 0.02 | 2.45 |
| CXCL13 | down | 0.04 | -2.50 |
| DROSHA | up | 0.04 | 1.32 |
| EP300 | up | 0.03 | 1.50 |
| EPHA2 | down | 0.01 | -11.42 |
| ERBB3 | down | 0.02 | -2.35 |
| GADD45B | down | 0.02 | -10.10 |
| HOXB2 | down | 0.03 | -4.26 |
| IKBKG | up | 0.03 | 1.71 |
| IL10 | down | 0.01 | -7.17 |
| IL1B | down | 0.03 | -10.47 |
| IL21 | down | 0.03 | -1.55 |
| IL4R | up | 0.03 | 3.12 |
| IL6 | down | 0.03 | -13.72 |
| IL8 | down | 0.02 | -30.61 |
| ING1 | down | 0.04 | -6.39 |
| JUN | down | 0.04 | -5.18 |
| KLF4 | down | 0.03 | -17.24 |
| LGALS3 | down | 0.03 | -9.39 |
| MAF | down | 0.03 | -1.64 |
| MAPK1 | up | 0.03 | 1.53 |
| MKI67 | down | 0.03 | -3.82 |
| MLL2 | up | 0.04 | 1.65 |
| MYC | down | 0.03 | -3.72 |
| NFKBIA | down | 0.04 | -2.06 |
| NUCB1 | up | 0.03 | 1.60 |
| PAX5 | up | 0.03 | 1.80 |
| PRDM1 | down | 0.03 | -7.93 |
| PRKCE | up | 0.03 | 2.36 |
| RAD51 | down | 0.03 | -7.24 |
| RUNX1 | up | 0.03 | 2.00 |
| S100A9 | down | 0.03 | -17.70 |
| SOCS2 | down | 0.04 | -2.86 |
| STAT4 | down | 0.03 | -8.22 |
| SYK | up | 0.03 | 1.65 |
| TNF | down | 0.03 | -3.41 |
| UBE2I | up | 0.03 | 1.46 |
| VEGFA | down | 0.03 | -11.82 |
| S. 7-7-14.5 | A TOTAL COLOR | 5177 | |